

Every minute matters when it comes to saving lives.

That's why companies are turning to quick-to-deploy life-saving equipment like the Evacusafe Excel Chair.

Evacusafe Excel chairs are designed specifically to help people with limited or impaired mobility.

Excel chairs are lightweight and easy to set up and use. With a friction-controlled track for slow and controlled descent, the Excel easily glides down stairs, fire escapes or flat surfaces, with up to 400 pounds on board, no lifting or manual handling required.

Once ground level has been reached, the Excel chair handles fold away and turn the device into a temporary wheelchair. With a supportive flat seat base, armrests and padded headrest the chair offers comfortable seating for extended periods of time.

When not in use, Evacusafe Excel chairs fold flat for hanging. Each chair ships with a set of wall hooks, an "Evacuation Chair" wall sign and a sturdy, long lasting protective cover.

Evacusafe Excel Chair Specifications:

Load Capacity	400 lbs.
Height (folded/open)	44.5"/52"
Depth (folded/open)	11.5"/39" (max)
Width	25.5"
Weight	30 lbs.
Length in use	56.25"

[Learn more.](#)
[Contact us today.](#)



Excel Evacuation Chair Benefits:

- Friction-controlled track helps the chair glide smoothly down stairs under single-person operation
- Four wheels make the chair easy to maneuver; also doubles as a temporary wheel chair
- Lightweight, easy to set up and use, providing swift and safe evacuation
- Supportive flat seat, cushioned head, arm and foot rests included at no additional cost.
- Foot operated brake and quick release safety harness for secure, safe operation