

# 180 AED cabinets

## Modern Metal Products

Modern Metal Products has been manufacturing cabinets in Owatonna Minnesota for over 85 years. We began our cabinet line with the fire extinguisher cabinet and have since started making many different sizes and styles of AED cabinets. MMP proudly stand by all of our quality manufactured cabinets.

## 180 AED Cabinet

The 180 series AED cabinets is a well designed versatile cabinet that is large enough for all AED's. Provides extra space for additional components and accessories.

- Made from 20 gage corrosive resistant steel. Stainless steel is available.
- Size is 17-1/2" x 17-1/2" x 7".
- Powder painted white to ensure a nice looking durable finish.
- Only requires one pull to open.
- Three trim styles to choose from:
  - Surfaces mount.
  - Recessed (ADA compliant)
  - Semi-recessed (ADA compliant)
- Choice of 3 features:
  - Cabinet only.
  - Cabinet with alarm.
  - Cabinet with alarm and strobe.
- Modern Metal Products AED cabinets only use MMP's 200900 alarm. This is the most versatile alarm in the indus-



Picue of model 180SM-1

**QUALITY CABINETS  
FROM QUALITY  
PEOPLE**

### **Modern Metal Products**

PO Box 247  
Owatonna, MN 55060  
Phone: 800-435-5544  
Fax: 507-451-0882  
E-mail: [gary@modern-metal.com](mailto:gary@modern-metal.com)  
[Www.modern-metal.com](http://www.modern-metal.com)

**MMP**  
MODERN METAL  
PRODUCTS