

EMERGENCY FIRST AID OXYGEN



FEATURES AND BENEFITS

- Meets FDA requirements – no prescription needed.
- Carried in durable shoulder bag – leaves both hands free
- Shoulder bag has a window to view contents gauge
- Permanently attached on/off handle – no wrench needed
- Pressure gauge – shows oxygen level in cylinder
- Refillable lightweight aluminum cylinder
- 4' kink proof hose

CONTENTS: 246-liter oxygen

USP FLOW: Regulator Flow = 1/2 LPM (LOW) to 15 LPM (High)

ENDURANCE: 246-liter (8 LPM) 30 minute supply **MATERIALS:** Aluminum tank; Chrome plated brass regulator

INCLUDES: 2 Masks with kink free hose

NOTE:

- 1) Review all product data sheets on proper use of equipment.
- 2) Oxygen Tank and Mask / Hose **MUST BE** pre-connected prior to putting equipment in service.
- 3) Most medical oxygen companies in your local area can refill your oxygen tank.