

Evacusafe Excel Evacuation Stair Chair

Features:

- 4 permanently deployed wheels with safety break on rear
- Friction controlled track
- Wide flat seat, fully open
- Weight bearing armrests
- Foot rest
- Padded head cushion with head restraint strap at no additional cost
- 3-point client support harness
- Anti-slip operator handles
- Integrated carry handles
- Additional leg strap
- Fold flat for wall mounting **Benefits:**
- 4 wheels aid stability and manoeuvrability on any landing or level ground
- Track allows for a safe, smooth, controlled desent down stair
- Open seat and weight bearing armrest make self-transfer for wheelchair users easier
- Flat seat provides familiar and comfy seating position
- Cushion offers support of head during the descent
- Over shoulder harness along with the arm rest helps those with weak core strength and offers additional safety over single belt
- Extra handling grip for operator during use
- Option to ascend stairs or over obstacles via integrated handles
- Extra safety strap to secure legs of persons with minimal to zero control
- Takes up minimal wall space

Whats Included:

- Protective dust cover
- Wall mounting hooks
- Wall sign
- Visual inspection sheet
- Operating guide
- Instructional DVD