


# NEW

## Major Bleeding Control Kit!

One kit with Surface Mount Cabinet \$149





### INCLUDES

- 1 Combat Tourniquet
- 1 Trauma Scissors
- 1 Bandage, 5x9
- 1 Trauma Dressing, 10x30
- 1 Kerlix
- 1 Roll of Medical Tape
- 1 Survival Blanket
- 8 Protective Rubber Gloves

Note: One time use kit, for major bleeding only.

### SURFACE MOUNTED CABINET

	INSIDE CABINET DIMENSIONS		
	Height	Width	Depth
<b>211SM</b>	11.25"	11.25"	3.25"



The 211SM is a surface mount, ADA compliant cabinet, measuring 11.25" (h) x 11.25" (w) x 3.25" (d). The unique mounting holes and free laying decal allow for a variety of door orientation options; left, right, top, or bottom. The compact cabinet provides endless storage and display possibilities and is ideal for small bleeding control or first aid kits.