


# **ZOLL AED 3**<sup>™</sup> The AED that leads the way

# **SMART**

# Real CPR Help®

- Real CPR Help guides rescuers in delivering high-quality CPR consistent with current guidelines
- Real-time feedback on depth and rate of compressions
- Full-color images and integrated voice, text, and visual prompts for support throughout the rescue
- Enhanced bar gauge allows rescuers to see when they are doing good CPR in real-time



# READY

# **Integrated Pediatric Rescue**

Universal CPR Uni-padz™ are always ready for both adult and child rescues

# **Low Cost of Ownership**

 Electrode pads and batteries last five years – fewer replacements and a lower total cost of ownership

# **WiFi Connectivity**

- Can automatically report AED readiness
- Quickly transfers case data

To learn more about the ZOLL AED 3<sup>™</sup>, call us at 1-(800) 804-4356 or visit us at www.zoll.com/aed3.

© 2017 ZOLL Medical Corporation. All rights reserved. Real CPR Help, Uni-padz, ZOLL AED 3, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries.

**ZOLL**。

Printed in USA. MCN PP 1710 0272