




☎ 844-277-2337

AFA-1 TLI

Surface Mount AED Cabinet



- The AED cabinet is a well designed versatile cabinet that is large enough for all AED's
- Provides extra space for additional components and accessories
- Made from 20 gage corrosive resistant steel stainless steel is available
- Size is 15" x 14" x 6 3/4
- Powder painted white to ensure a nice looking durable finish
- Only requires one pull to open
- Surfaces mount

INCLUDED

- Door Alarm (9-volt)
- Includes Steel V-Shaped Wall Sign
- AED Onsite Window Decal
- How to CPR/AED Poster

Sudden Cardiac Arrest Treatment

CPR / AED Guidelines

If the Victim displays the following signs:

- Unconscious
- Not Responsive
- Not breathing normally

Steps you should do:

- Call 911 and get help
- Obtain an AED
- Start hands-only compression CPR - place hands in center of the chest and press hard (2" deep) and fast (about 100 compressions per minute)
- When an AED arrives, turn it on and follow all voice prompts until 911 arrives

