

LIFEPAK CR2**TOUGHcase**

Designed specifically for the LIFEPAK CR2, the TOUGHcase™ hard case provides rugged protection while maintaining rapid access to critical device functions during storage, transport and emergency response. Engineered for demanding indoor and outdoor environments, the case was validated through environmental, usability and performance testing to support real-world responder workflows.

**Built to carry the LIFEPAK TOUGH legacy**

- Minimal storage volume that profiles a CR2 shape
- Ensures uninterrupted Wi-Fi and cellular connectivity of the CR2
- User can confirm CR2 readiness from case exterior without opening the case
- Case can withstand heavy force without crushing or breaking and can support 300 lbs (136 kg) of weight on it.
- Can be opened via one latch with one hand
- Easy to clean and can be wiped clean by hand
- IP Rating of 68
- Waterproof to a depth of 2 meters for up to 60 minutes with a CR2 inside
- The case is buoyant with up to 10 lbs. (4.6 kg) of contents, including the CR2
- Case only weighs 4 lbs. (1.8 kg)

AED users should be trained in CPR and in the use of the AED. AEDs which are for public use/trained users can also be used by professional users.

The information presented is intended to demonstrate Stryker's product offerings. Refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use, or have applied for the following trademarks or service marks: LIFEPAK, LIFEPAK Tough, TOUGHcase and Stryker. All other trademarks are trademarks of their respective owners or holders.